



Hefei Sada Medical Equipment Co., Ltd.

Address: No. 1412, JingDian Plaza, Economic Development Zone, Hefei, Anhui, China

Email: kevin@sadamedical.com, Website: www.sadamedical.com

Tel/Fax: +86-551-63879990, Mobile / WhatsApp /: +86-15955157922

Portable Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure; Hand wheel sliding door 2. Digital display of working status, touch key for setting Temp. and time 3. Electric heated 4. Double scale indication pressure gauge 5. Automatically shut off with beep reminding after sterilization <p>Technical Data:</p> <p>Working Pressure: 0.14-0.165Mpa Working Temperature: 105°C-126°C Timer: 0-99 minutes Power: AC220V, 50 HZ, 2kw</p>	PA-ND18	18L Ø280*H260mm 1pc Sterilizing Basket	400*400*610	12
		PA-ND24	24L Ø280*H390mm 1pc Sterilizing Basket	400*400*730	14
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure; Hand wheel sliding door 2. Indicator lights indicate working status 3. Electric heated 4. Double scale indication pressure gauge 5. Temp. and time can be set by knobs 6. Automatically shut off with beep reminding after sterilization <p>Technical Data:</p> <p>Working Pressure: 0.14-0.165Mpa Working Temperature: 105°C-126°C Timer: 0-60 minutes Power: AC220V, 50 HZ, 2KW</p>	PA-NJ18	18L Ø280*H260mm 1pc Sterilizing Basket	400*400*610	12
		PA-NJ24	24L Ø280*H390mm 1pc Sterilizing Basket	400*400*730	14

Please find your preferred models and send email to us for request Latest Best Price



Hefei Sada Medical Equipment Co., Ltd.

Address: No. 1412, JingDian Plaza, Economic Development Zone, Hefei, Anhui, China

Email: kevin@sadamedical.com, Website: www.sadamedical.com

Tel/Fax: +86-551-63879990, Mobile / WhatsApp /: +86-15955157922

Portable Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure; Hand wheel sliding door 2. Electric or gas heated 3. Double scale indication pressure gauge 4. Easy to operate, safe and reliable <p>Technical Data:</p> <p>Working Pressure: 0.14-0.165Mpa Working Temperature: 105°C-126°C Power: 220V, 50 HZ, 2KW</p>	PA-NM12	12L Ø240*H200mm 1pc Sterilizing Basket	340*340*500	10
		PA-NM18	18L Ø280*H265mm 1pc Sterilizing Basket	400*400*530	11
		PA-NM24	24L Ø280*H390mm 1pc Sterilizing Basket	400*400*630	13
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Electric or gas heated 3. Double scale indication pressure gauge 4. Easy to operate, safe and reliable <p>Technical Data:</p> <p>Working Pressure: 0.14-0.165Mpa Working Temperature: 105°C-126°C Power: 220V, 50 HZ, 2KW</p>	PA-BM12	12L Ø240*H200mm 1pc Sterilizing Basket	340*340*400	11
		PA-BM18	18L Ø280*H260mm 1pc Sterilizing Basket	390*390*410	14
		PA-BM24	24L Ø280*H390mm 1pc Sterilizing Basket	390*390*530	15

Please find your preferred models and send email to us for request Latest Best Price



Hefei Sada Medical Equipment Co., Ltd.

Address: No. 1412, JingDian Plaza, Economic Development Zone, Hefei, Anhui, China

Email: kevin@sadamedical.com, Website: www.sadamedical.com

Tel/Fax: +86-551-63879990, Mobile / WhatsApp /: +86-15955157922

Portable Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Electric or gas heated 3. Double scale indication pressure gauge 4. Easy to operate, safe and reliable <p>Technical Data: Working Pressure: 0.145-0.165Mpa Working Temperature: 126°C-129°C Power: 220V, 50 HZ, 2KW</p>	YX-280B-18L	18L Ø280*H243mm 1pc Sterilizing Basket	380*380*400	15
		YX-280B-24L	24L Ø280*H383mm 1pc Sterilizing Basket	380*380*510	17
		YX-280B-30L	30L Ø280*H423mm 1pc Sterilizing Basket	380*380*620	19
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Electric heated 3. Double scale indication pressure gauge 4. With on/off switch & drain water tap, easy to operate <p>Technical Data: Working Pressure: 0.145-0.165Mpa Working Temperature: 126°C-129°C Power: 220V, 50 HZ, 2KW-2.5KW</p>	YX-280A-18L	18L Ø280*H243mm 1pc Sterilizing Basket	380*380*510	16
		YX-280A-24L	24L Ø280*H383mm 1pc Sterilizing Baskets	380*380*620	17
		YX-280A30	30L Ø280*H423mm 1pc Sterilizing Basket	380*380*680	19

Please find your preferred models and send email to us for request Latest Best Price



Hefei Sada Medical Equipment Co., Ltd.

Address: No. 1412, JingDian Plaza, Economic Development Zone, Hefei, Anhui, China

Email: kevin@sadamedical.com, Website: www.sadamedical.com

Tel/Fax: +86-551-63879990, Mobile / WhatsApp /: +86-15955157922

Portable Autoclave Sterilizer Catalog

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Electric heated 3. Double scale indication pressure gauge 4. Temp. and time can be set by knobs 5. With on/off switch 6. Indicator lights indicate working status <p>Technical Data:</p> <p>Working Pressure: 0.145-0.165Mpa Working Temperature: 126°C-129°C Timer: 0-80 minutes Power: 220V, 50 HZ, 2KW-2.5KW</p>	YX-280D-18L	18L Ø280*H243mm 1pc Sterilizing Basket	380*380*510	16
		YX-280D-24L	24L Ø280*H383mm 1pc Sterilizing Basket	380*380*620	17
		YX-280D-30L	30L Ø280*H423mm 1pc Sterilizing Basket	380*380*680	19
	<p>Features:</p> <ol style="list-style-type: none"> 1. Full stainless steel structure 2. Dual scale numerical pressure gauge indicate temperature and pressure 3. Digital display of working status, touch key for setting time and temperature 4. Power emergency control switch, Over heat protection 5. Automatic cut off power if the water is not enough <p>Technical Data:</p> <p>Working Pressure: 0.145-0.165Mpa Working Temperature: 126°C-129°C Timer: 0-999 minutes Power: 220V, 50 HZ, 2KW</p>	YX-280DG-18L	18L Ø280*H243mm 1pc Sterilizing Basket	380*380*510	16
		YX-280DG-24L	24L Ø280*H383mm 1pc Sterilizing Basket	380*380*620	17
		YX-280DG-30L	30L Ø280*H423mm 1pc Sterilizing Basket	380*380*680	19

Please find your preferred models and send email to us for request Latest Best Price