

Table Top Sterilizer Autoclave Catalogue



Picture	Description	Model No.	Volume	Product Dim. (mm)	Package Size (mm)
	<p>Features:</p> <ol style="list-style-type: none"> Indicator light indicates working state. For 4~6 minutes rapidly sterilizing. Sterilizing temperature and time can be preset. Steam-water inner circulation system:no steam discharge,and the environment for sterilizing will be clean and dry. Over temperature & over pressure auto-protection. Safe protection of water lacking. Exhaust the cool air automatically. With three stainless steel sterilizing plates. The chamber of the sterilizer is made of stainless steel. Automatically shut off with beep reminding after sterilization. <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~134°C Timer: 0~60min Power: 1.5KW / AC220V 50HZ</p>	TS-AJ20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x384	700*580*500
		TS-AJ24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	800*580*500
	<p>Features:</p> <ol style="list-style-type: none"> Indicator light indicates working state. For 4~6 minutes rapidly sterilizing. Sterilizing temperature and time can be preset. Steam-water inner circulation system:no steam discharge,and the environment for sterilizing will be clean and dry. Over temperature & over pressure auto-protection. Safe protection of water lacking. Exhaust the cool air automatically. Door safety lock system. With four stainless steel sterilizing plates. The chamber of the sterilizer is made of stainless steel. Automatically shut off with beep reminding after sterilization. <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~134°C Timer: 0~60min Power: 2.7KW / AC220V 50HZ</p>	TS-C35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	820*680*650
		TS-C50	50L Ø340*H550mm 4pcs sterilizing plate	770x600x540	880*730*700

Table Top Sterilizer Autoclave Catalogue



Picture	Description	Model No.	Volume	Product Dim. (mm)	Package Size (mm)
	<p>Features:</p> <ol style="list-style-type: none"> 1. For 4~6 minutes rapidly sterilizing. 2. Digital display of working status,touch type key. 3. With 3 fixed cycles of sterilizing program and adjustable program for user. 4. The whole process of water adding,temperature rising,sterilizing,drying steam discharging controlled automatically. 5. Steam-water inner circulation system:no steam discharge,and the environment for sterilizing will be clean and dry. 6. Exhaust the cool air automatically. 7. Safe protection of water lacking. 8. Door safety lock system. 9. With three stainless steel sterilizing plates. 10. The chamber of the sterilizer is made of stainless steel. 11. Automatically shut off with beep reminding after sterilization. 12. With drying function. <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~99min Power: 1.5-2.7KW / AC220V 50HZ</p>	TS-CD20	20L Ø250*H420mm 3pcs sterilizing plate	480x480x380	700*580*500
		TS-CD24	24L Ø250*H520mm 3pcs sterilizing plate	580x480x384	800*580*500
		TS-CD35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	810*680*650
		TS-CD50	50L Ø340*H550mm 4pcs sterilizing plate	770x600x540	800*730*700

Table Top Sterilizer Autoclave Catalogue




Picture	Description	Model No.	Volume	Product Dim. (mm)	Package Size (mm)
	<p>Features:</p> <ol style="list-style-type: none"> Standard Class B with three times vacuum and drying, the remaining temperature of sterilized instrument is less than 0.2%. LCD display working status, touch keys. The vacuum data can reach -0.8bar, it fits various instrument, including package, unpackaged, solid, hollow, multi-apertures and inset pipes. It makes sure the instrument can be sterilized completely. Fully computer controlled LCD display bright show. The interface is modularized panel for easy operation. It is installed with BOWIE & DICK which measures the penetration of water steam. It has vacuum text which can make sure the ability of texting vacuuming. Fast and independent steam generator. Built-in mini printer which can record the sterilizer information. Safety valve, Safety lock door system, Pressure or temperature overload protection, Alarm system. Can inactivate HIV, HBV, mad cow virus and bacillus. <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~134°C Timer: 0~99min Power: 2.0-2.7KW / AC220V 50HZ</p>	TS-DDV18	18L Ø250*H360mm 3pcs sterilizing plate	560x470x400	655*560*495
		TS-DDV24	24L Ø250*H470mm 3pcs sterilizing plate	675x470x400	765*560*495
		TS-DDV35	35L Ø300*H500mm 4pcs sterilizing plate	730x550x510	810*680*650
		TS-DDV50	50L Ø340*H550mm 4pcs sterilizing plate	770x600x540	880*730*700

Table Top Sterilizer Autoclave Catalogue

Picture	Description	Model No.	Volume	Product Dim. (mm)	Package Size (mm)
	<p>Features:</p> <ol style="list-style-type: none"> 1. Sterilizing course: Automatic sterilization controlled by computer, easy to operate. 2. Maximum temperature: up to 134°C, suit for 4-6 minutes rapidly sterilizing. 3. Safe protection of water lacking. 4. With over-temperature, over-pressure auto-protect device 5. After sterilization it will automatically cut the power and then alarm 6. Fully stainless steel structure. <p>Technical Data: Working pressure: 0.22Mpa Working temperature: 134°C Adjustment of temperature: 105~135°C Timer: 0~60min Power: 1.5KW / AC220V 50HZ</p>	<p>TS-B24</p>	<p>24L Ø270*H410mm 3pcs sterilizing plate</p>	<p>560x350x500</p>	<p>650*410*570</p>