




## Automatic Vertical Pressure Steam Autoclave

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. Fully stainless steel structure.</li> <li>2. Hand wheel type of quick-open door structure.</li> <li>3. Door safety lock system.</li> <li>4. Automatically controlled by computer.</li> <li>5. Digital display of working status, touch type key.</li> <li>6. Over temperature and over pressure auto-protection.</li> <li>7. Steam- water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.</li> <li>8. Safe protection of water lacking.</li> <li>9. Self-inflating type seal.</li> <li>10. Automatically shut off with beep reminding after sterilization.</li> <li>11. With two stainless steel sterilizing baskets.</li> <li>12. With drying system. Printer function (USD340.00 if need).</li> </ol> <p><b>Technical Data:</b></p> <p>Working pressure: 0.22Mpa            Working temperature: 134°C            Adjustment of temperature: 105~134°C            Timer: Timer: 0~99 min or 0-99 hour 59 min            Power:2.5-4.5KW / AC 220V±22V, 50Hz±1Hz</p>	VA-SA35	35L Ø318*H450mm 2pcs sterilizing baskets	600*410*1030	77
		VA-SA50	50L Ø340*H550mm 2pcs sterilizing baskets	650*410*1140	83
		VA-SA75	75L Ø400*H600mm 2pcs sterilizing basket	690*470*1140	98
		VA-SA100	100L Ø440*H650mm 2pcs sterilizing baskets	730*510*1270	110


## Automatic Vertical Pressure Steam Autoclave

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. Fully stainless steel structure.</li> <li>2. Hand wheel type of quick-open door structure.</li> <li>3. Door safety lock system.</li> <li>4. Digital display of working status, touch type key.</li> <li>5. Auto discharge the cool air, and steam discharging automatically after sterilization</li> <li>6. Over temperature and over pressure auto-protection.</li> <li>7. Safe protection of water lacking.</li> <li>8. Self-inflating type seal.</li> <li>9. Automatically shut off with beep reminding after sterilization.</li> <li>10. With two stainless steel sterilizing baskets.</li> <li>11. Drying system.(USD170 if need)</li> </ol> <p><b>Technical Data:</b></p> <p>Working pressure: 0.22Mpa            Working temperature: 134℃            Adjustment of temperature: 105~134℃            Timer: Timer: 0~99 min or 0-99 hour 59 min            Power: 2.5-4.5KW / AC220V 50HZ</p>	VA-SD35	35L Ø318*H450mm 2pcs sterilizing baskets	450*450*1010	56
		VA-SD50	50L Ø340*H550mm 2pcs sterilizing baskets	510*470*1130	68
		VA-SD75	75L Ø400*H600mm 2pcs sterilizing baskets	560*560*1120	80
		VA-SD100	100L Ø440*H650mm 2pcs sterilizing baskets	540*560*1250	85


## Automatic Vertical Pressure Steam Autoclave

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. Automatically controlled by computer, digital display of working status, touch type key;</li> <li>2. Dynamically showed the temperature, pressure time and vacuum rate, and the full cycle;</li> <li>3. Fully automatically route for water inlet, pulse vacuum, heating, sterilizing, draining, exhausting, vacuum drying;</li> <li>4. 3 times pulse vacuum, which make the vacuum rate as -0.081 Mpa, and ensure the high extracting of inner air from dressing, A, B type hollowed appliance;</li> <li>5. With 4 fixed program for appliance, medical dressing, rubbers and fluid, and the user can also set the program as they need;</li> <li>6. Equipped with the PT/TT testing connector, easy to run the testing;</li> <li>7. Multiple pulse vacuum drying after sterilizing, ensure the nice sterilizing;</li> <li>8. Equipped with 0.22 μm bacteria filter, ensure the sterile air enter the inner vacuum environment, avoid the second contamination;</li> <li>9. Hand wheel type of quick-open door structure, self-inflating type seal;</li> <li>10. Adopted with immersed heating, with high heating efficiency;</li> <li>11. Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry;</li> <li>12. Automatically shut off with beep reminding after sterilization;</li> <li>13. Fully high quality stainless, equipped with double stainless steel baskets;</li> <li>14. Printer function (USD340.00 if need).</li> </ol> <p><b>Technical Data:</b>            Working pressure: 0.22Mpa            Working temperature: 134℃            Adjustment of temperature: 0~134℃            Timer: 1~99min            Power: 1.5-4.5KW / AC220V 50HZ</p>	VA-SV50	50L Ø340*H550mm 2pcs sterilizing baskets	690*470*1140	105
		VA-SV75	75L Ø400*H600mm 2pcs sterilizing baskets	730*510*1270	125
		VA-SV100	100L Ø440*H650mm 2pcs sterilizing baskets	770*550*1350	140

## Automatic Vertical Pressure Steam Autoclave

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
 <p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. All material is SUS304 stainless steel.</li> <li>2. Digital display of working status, touch type key.</li> <li>3. Auto discharge the cool air, and steam discharging automatically after sterilization.</li> <li>4. Automatically shut off with beep reminding after sterilization.</li> <li>5. With two stainless steel baskets.</li> </ol> <p><b>Technical Data:</b></p> <p>Working pressure: 0.22Mpa            Working temperature: 134℃            Adjustment of temperature: 105~134℃            Timer: 0~99 min or 0-99 hour 59 min            Power: 2.5-6KW / AC220V 50HZ</p>	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>1. All material is SUS304 stainless steel.</li> <li>2. Digital display of working status, touch type key.</li> <li>3. Auto discharge the cool air, and steam discharging automatically after sterilization.</li> <li>4. Automatically shut off with beep reminding after sterilization.</li> <li>5. With two stainless steel baskets.</li> </ol> <p><b>Technical Data:</b></p> <p>Working pressure: 0.22Mpa            Working temperature: 134℃            Adjustment of temperature: 105~134℃            Timer: 0~99 min or 0-99 hour 59 min            Power: 2.5-6KW / AC220V 50HZ</p>	VA-FD35	35L Ø318*H450mm 2pcs sterilizing baskets	480*460*850	42
		VA-FD50	50L Ø340*H550mm 2pcs sterilizing baskets	520*520*980	50
		VA-FD75	75L Ø400*H600mm 2pcs sterilizing baskets	560*560*980	70
		VA-FD100	100L Ø440*H650mm 2pcs sterilizing baskets	590*590*1080	85
		VA-FD120	120L Ø480*H660mm 2pcs sterilizing baskets	600*640*1140	100
		VA-FD150	150L Ø510*H740mm 2pcs sterilizing baskets	670*690*1140	110

## Automatic Vertical Pressure Steam Autoclave

Picture	Description	Model No.	Volume	Product Dim. (mm)	Net Weight (kg)
	<p><b>Features:</b></p> <ol style="list-style-type: none"> <li>7" capacitive touch screen, fully automatically controlled, easy to operate on.</li> <li>Elevating type door, which is more convenient to place and fetch the objects from the autoclave.</li> <li>The whole procedure from water inlet, pulse vacuum, heat, sterilizing, water outlet, exhausting to vacuum drying, which are all automatically by one-key operation.</li> <li>Pulse vacuum, the limit is-0.080Mpa, which ensure the sterilizing more completely</li> <li>Separated steam generator, which is more energy saving by comparing to the normal heating and prevent the sterilizing objects from been oxidized.</li> <li>There are several programs loaded in the autoclaves for different requirement of different objects</li> <li>With PT-TT testing port, easy to run the testing</li> <li>Equipped with Germany SARTORIUS Sterile ventilation filter MIDISART 2000 0.2 um, which ensure the air is sterile, prevent the contamination.</li> <li>Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.</li> </ol> <p><b>Technical Data:</b>            Working pressure: 0.22Mpa            Working temperature: 134℃            Adjustment of temperature: 105~134℃            Timer: 1~99min            Power: 2.4-3.6KW / AC220V 50HZ</p>	VA-SVN60	60L Ø396×490mm	686*556*913	180
		VA-SVN80	80L Ø396×650mm	686*556*1035	200
		VA-SVN100	100L 504×404×502mm	840*600*970	270
		VA-SVN120	120L 504×404×590mm	840*600*1060	295