



## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Self Sealing Sterilization Pouch (200pcs/Box)

Item No.	Size MM x MM	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
			L	W	H				
SSP.057-130	57x130	60box/ctn	0.37	0.3	0.3	11.5	11	Disinfecting Type: Steam and EO Shelf Life:3 years Material: Medical grade paper + Film Packing: 200pcs/Box Quality Certification: TUV Payment:TT MOQ: 100000PCS Sample Free BRAND CAN BE CUSTOMIZED	
SSP.070-260	70x260	20box/ctn	0.31	0.3	0.3	9.5	9		
SSP.090-165	90x165	20box/ctn	0.39	0.2	0.3	8	7.5		
SSP.090-260	90x260	20box/ctn	0.39	0.3	0.3	12	11.5		
SSP.135-260	135x260	10box/ctn	0.29	0.3	0.3	9	8.5		
SSP.135-280	135x280	10box/ctn	0.31	0.3	0.3	10	9.5		
SSP.190-360	190x360	5box/ctn	0.33	0.4	0.2	9	8.5		
SSP.250-370	250x370	5box/ctn	0.4	0.3	0.27	12	11.5		
SSP.260-410	260x410	5box/ctn	0.43	0.3	0.3	13.5	13		
SSP.305-410	305x410	5box/ctn	0.43	0.3	0.3	16	15.5		
SSP.305-430	305x430	5box/ctn	0.45	0.3	0.3	16.5	16		
SSP.300-490	300x490	5box/ctn	0.51	0.3	0.3	18.5	18		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Sterilization Flat Roll (200meter/Roll)

Item No.	Size MM x M	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SSR.050-200	50x200	12roll/ctn	0.42	0.2	0.3	14.5	14	Disinfecting Type: ETO/STEAM/GAMMA RAY Sterilization Shelf Life:3 years Material: Medical grade paper + Film Packing: 200 Meters/Roll Quality Certification: ISO/CE Payment:TT MOQ: 200 Rolls Sample Price: USD\$10/Roll OEM Order Acceptable	
SSR.055-200	55x200	12roll/ctn	0.42	0.2	0.4	15.5	15		
SSR.075-200	75x200	8roll/ctn	0.42	0.2	0.3	14.5	14		
SSR.100-200	100x200	6roll/ctn	0.42	0.2	0.3	14.5	14		
SSR.150-200	150x200	4roll/ctn	0.42	0.2	0.3	14.5	14		
SSR.200-200	200x200	2roll/ctn	0.42	0.2	0.2	9.5	9		
SSR.250-200	250x200	2roll/ctn	0.42	0.2	0.3	12	11.5		
SSR.300-200	300x200	2roll/ctn	0.42	0.2	0.3	14.5	14		
SSR.350-200	350x200	1roll/ctn	0.22	0.2	0.4	8.5	8		
SSR.400-200	400x200	1roll/ctn	0.22	0.2	0.4	9.5	9		
SSR.450-200	450x200	1roll/ctn	0.22	0.2	0.5	11	10.5		
SSR.500-200	500x200	1roll/ctn	0.22	0.2	0.5	12	11.5		
SSR.600-200	600x200	1roll/ctn	0.22	0.2	0.6	14.5	14		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Sterilization Gusseted Roll (100meter/Roll)

Item No.	Size MM x M	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SGR.075-100	75x100	8rolls/tn	0.4	0.2	0.3	10	9.5	Disinfecting Type: ETO/STEAM/GAMMA RAY Sterilization Shelf Life:3 years Material: Medical grade paper + Film Packing: 200 Meters/Roll Quality Certification: ISO/CE Payment:TT MOQ: 200 Rolls Sample Price: USD\$10/Roll OEM Order Acceptable	
SGR.100-100	100x100	6rolls/tn	0.4	0.2	0.3	10	9.5		
SGR.150-100	150x100	4rolls/tn	0.4	0.2	0.3	10	9.5		
SGR.200-100	200x100	2rolls/tn	0.4	0.2	0.3	7.5	7		
SGR.250-100	250x100	2rolls/tn	0.4	0.2	0.3	8	7.5		
SGR.300-100	300x100	2rolls/tn	0.4	0.2	0.3	10	9.5		
SGR.350-100	350x100	2rolls/tn	0.4	0.2	0.4	11	10.5		
SGR.400-100	400x100	2rolls/tn	0.4	0.2	0.4	12	11		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### TYVEK Sterilization Roll (70meters/Roll)

Item No.	Size MM x M	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
			L	x W	x H				
TVKSSR.075-70	75 x70	16rolls/ctn	0.34	0.3	0.3	10.5	10	Disinfecting Type: ETO/STEAM/GAMMA RAY Sterilization Shelf Life:3 years Material: Medical grade paper + Film Packing: 200 Meters/Roll Quality Certification: ISO/CE Payment:TT MOQ: 100 Rolls Sample Price: USD\$10/Roll OEM Order Acceptable	
TVKSSR.100-70	100 x70	12rolls/ctn	0.34	0.3	0.3	10.5	10		
TVKSSR.150-70	150 x70	8rolls/ctn	0.34	0.3	0.3	10.5	10		
TVKSSR.200-70	200 x70	4rolls/ctn	0.34	0.3	0.3	7.5	7		
TVKSSR.250-70	250 x 70	4rolls/ctn	0.34	0.3	0.3	9	8.5		
TVKSSR.300-70	300x70	4rolls/ctn	0.34	0.3	0.3	10.5	10		
TVKSSR.350-70	350x70	4rolls/ctn	0.34	0.3	0.4	12.5	12		
TVKSSR.400-70	400x70	4rolls/ctn	0.34	0.3	0.4	14	13.5		
TVKSSR.450-70	450x70	4rolls/ctn	0.34	0.3	0.5	12.5	12		
TVKSSR.500-70	500x70	2rolls/ctn	0.34	0.3	0.5	9	8.5		
TVKSSR.600-70	600x70	2rolls/ctn	0.34	0.3	0.6	10.5	10		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### 60gsm Non-woven Wrap

Item No.	Size CM x CM	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
			L	W	H				
SSNW3030	30x30	1000pc/ctn	0.62	0.3	0.2	6	5.5	Disinfecting Type: ETO/STEAM/GAMMA RAY Sterilization Shelf Life:3 Years Material: 100%Polypropylene Thickness: 0.01mm-1mm Weight: 45GSM/60GSM Packing: 100 PCS/Bag Quality Certification: ISO13458/CE Payment:TT MOQ: 1000PCS Sample Price: USD\$10/Roll OEM Order Acceptable	
SSNW4040	40x40	500pc/ctn	0.42	0.4	0.2	5.5	5		
SSNW4545	45x45	500pc/ctn	0.47	0.5	0.2	6.5	6		
SSNW5050	50x50	500pc/ctn	0.52	0.5	0.2	8	7.5		
SSNW6060	60x60	250pc/ctn	0.62	0.3	0.2	6	5.5		
SSNW7070	70x70	250pc/ctn	0.72	0.4	0.2	8	7.5		
SSNW7575	75x75	250pc/ctn	0.77	0.4	0.2	9	8.5		
SSNW8080	80x80	150pc/ctn	0.82	0.4	0.2	6.5	6		
SSNW9090	90x90	150pc/ctn	0.92	0.3	0.2	8	7.5		
SSNW100100	100x100	150pc/ctn	1.02	0.4	0.2	9.5	9		
SSNW120120	120x120	100pc/ctn	1.22	0.3	0.2	9	8.5		



## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### 60GSM CREPE PAPER

Item No.	Size CM x CM	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
			L	W	H				
SWP3030	30x30	2000pc/ctn	0.62	0.3	0.2	9	8.5	Disinfecting Type: ETO, Steam, Gamma Irradiation or E-Beam Sterilization Methodology Shelf Life:3 Years Material: 100% Wood Pulp, NOT Petroleum Coating: Uncoated, Water Proof Color: Blue/Green/White Packing: 500 Sheets/Carton Quality Certification: CE Payment:TT MOQ: 1000PCS Sample Price: FREE OEM Order Acceptable	
SWP4040	40x40	1000pc/ctn	0.42	0.4	0.2	8	7.5		
SWP4545	45x45	1000pc/ctn	0.47	0.5	0.2	10	9.5		
SWP5050	50x50	1000pc/ctn	0.52	0.5	0.2	12	11.5		
SWP6060	60x60	500pc/ctn	0.62	0.3	0.2	9	8.5		
SWP7070	70x70	500pc/ctn	0.72	0.4	0.2	12	11.5		
SWP7575	75x75	500pc/ctn	0.77	0.4	0.2	13.5	13		
SWP8080	80x80	300pc/ctn	0.82	0.4	0.1	9.5	9		
SWP9090	90x90	300pc/ctn	0.92	0.3	0.2	12	11.5		
SWP100100	100x100	200pc/ctn	1.02	0.4	0.1	10	11.5		
SWP120120	120x120	100pc/ctn	1.22	0.3	0.1	7.5	7		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China


### DENTAL BIBS

Item No.	Size CM x CM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SDB	33 x 45	500pcs/ctn	0.34	0.3	0.3	3.5	3	Shelf Life:3 Years Material: 2ply Paper+1ply PE Film Packing: PE Bag Quality Certification: CE Payment:TT MOQ: 10000PCS Sample Price: FREE OEM Order Acceptable	

### Steam Indicator Tape

Item No.	Size Width(MM) x Length(M)	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SSIT1250	12x50	180rolls/tn	0.42	0.4	0.3	14.5	14	Disinfecting Type: Steam Sterilization Methodology Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: CE Payment:TT MOQ: 500PCS Sample Price: USD\$10/Roll OEM Order Acceptable	
SSIT1950	19x50	117rolls/tn	0.42	0.4	0.3	14.5	14		
SSIT2550	25x50	90rolls/tn	0.42	0.4	0.3	14.5	14		

### EO Indicator Tape

Item No.	Size Width(MM) x Length(M)	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SEIT1250	12x50	180rolls/tn	0.42	0.4	0.3	14.5	14	Disinfecting Type: EO Sterilization Methodology Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: CE Payment:TT MOQ: 500PCS Sample Price: USD\$10/Roll OEM Order Acceptable	
SEIT1950	19x50	117rolls/tn	0.42	0.4	0.3	14.5	14		
SEIT2550	25x50	90rolls/tn	0.42	0.4	0.3	14.5	14		






## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Plasma Indicator Tape

Item No.	Size Width(MM) x Length(M)	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SPIT1250	12x50	95rolls/tn	0.35	0.3	0.3	6.5	6	Disinfecting Type: Plasma Sterilization Methodology Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: CE Payment:TT MOQ: 500PCS Sample Price: USD\$10/Roll OEM Order Acceptable	
SPIT1950	19x50	60rolls/tn	0.35	0.3	0.3	6.5	6		
SPIT2550	25x50	45rolls/tn	0.35	0.3	0.3	6.4	6		




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Steam Indicator Strip/Class 4

Item No.	Size MM x MM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SAIS20-140	20x140mm 20x70cc/pc, two pc/one unit	250 unit/box; 50boxes/ctn	0.5	0.2	0.2	9	8.5	Usage: Test Sterilization Procedure Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: ISO13485/CE MOQ: 10000PCS Sample Price: FREE	




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Steam Indicator Strip/Class 5

Item No.	Size MM x MM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SAIS20-140	20x140mm 20x70cc/pc, two pc/one unit	250 unit/box; 50boxes/ctn	0.5	0.2	0.2	9	8.5	Usage: Test Sterilization Procedure Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: ISO13485/CE MOQ: 10000PCS Sample Price: FREE	



## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Steam Indicator Strip/Class 6

Item No.	Size MM x MM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SAIS20-140	20x140mm 20x70cc/pc, two pc/one unit	250 unit/box; 50boxes/ctn	0.5	0.2	0.2	9	8.5	Usage: Test Sterilization Procedure Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: ISO13485/CE MOQ: 10000PCS Sample Price: FREE	




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519


Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China


### EO Indicator strip

Item No.	Size MM x MM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SAEI15-100	15x100mm	250 PCS/box; 100boxes/ctn	0.45	0.2	0.2	8.5	8	Usage: Test Sterilization Procedure Shelf Life:3 Years Material: Medical Grade Paper Quality Certification: ISO13485/CE MOQ: 10000PCS Sample Price: FREE	

### BD Test Pack

Item No.	Size MM x MM	Packing	CTN(meter) L x W x H			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SBDT115125	115*125mm	100PCS/box;	0.52	0.3	0.2	10	9.5	Usage: Test Sterilization Procedure Shelf Life:3 Years Material: Special BD Test Paper Quality Certification: CE Packing: 1PC/PE Pouch, 100PCS / CTN MOQ: 1000PCS Color Change: From Yellow to Black Sample Price: USD\$10	

### Pressure Steam Sterilization Biological Indicator

Item No.	Size	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SSBIN	1 x 4.6cm	25PCS/box; 20 boxes/CTN	0.32	0.3	0.4	6.28	6	Disinfecting Type: Steam Sterilizer Shelf Life:2 Years Material: Plastic Bottle Quality Certification: CE Packing: 25PCS/INNER BOX MOQ: 20 INNER BOXES Sample Price: USD\$45	




## Hefei Sada Medical Equipment Co., Ltd.

Tel/Fax: +86-551-63879990, Mobile / WhatsApp / WeChat : +86-18056076519

Email: vivian@sadamedical.com, Website: <https://www.sadamedical.com>

No. 12629, Jingdian Plaza, Fanhua Avenue, Economic Development Zone, Hefei, Anhui, China

### Hydrogen Peroxide Plasma Biological Indicator

Item No.	Size	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SPBIN	1 x 4.6cm	25PCS/box; 20 boxes/CTN	0.32	0.3	0.4	6.28	6	Disinfecting Type: Plasma Sterilizer Shelf Life: 2 Years Material: Plastic Bottle Quality Certification: CE Packing: 25PCS/INNER BOX MOQ: 20 INNER BOXES Sample Price: USD\$45	

### EO Sterilizer Biological Indicator

Item No.	Size	Packing	CTN(meter)			G.W.(KG)	N.W.(KG)	FEATURES	PHOTO
SEOBIN	1 x 4.6cm	25PCS/box; 20 boxes/CTN	0.32	0.3	0.4	6.28	6	Disinfecting Type: EO Sterilizer Shelf Life: 2 Years Material: Plastic Bottle Quality Certification: CE Packing: 25PCS/INNER BOX MOQ: 20 INNER BOXES Sample Price: USD\$45	